

SPECIFICATIONS




RAPID EXCHANGE PTCA CATHETER

- Usable length: 140cm
- Proximal shaft diameter: 1.7F (0.58mm)
- Distal shaft diameter: 2.3F (0.78mm) for Ø1.5 & 2.0mm
2.6F (0.86mm) for Ø2.5 & 3.5mm
2.8F (0.95mm) for Ø4.0 mm
- Minimum guiding catheter diameter: 5F (0.058")
- Maximum guide diameter: 0.014" (0.356mm)
- Radio-opaque markers: 2 of 0.75mm

	Pressure Bars	Balloon size (mm) (at 37°C)					
		Ø1.5 mm	Ø2.0 mm	Ø2.5 mm	Ø3.0 mm	Ø3.5 mm	Ø4.0 mm
Nominal pressure ✓	2	1.39	1.85	2.32	2.78	3.24	3.71
	4	1.44	1.92	2.40	2.88	3.36	3.84
	6	1.50	2.00	2.50	3.00	3.50	4.00
	8	1.54	2.05	2.56	3.07	3.58	4.09
	10	1.57	2.09	2.61	3.13	3.65	4.17
	12	1.59	2.12	2.65	3.18	3.71	4.24
Rated Burst Pressure* (RBP)	14	1.62	2.16	2.70	3.24	3.78	4.32
	16	1.65	2.20	2.75	3.30	3.85	4.40
	18	1.68	2.24	2.80	3.36	3.92	4.48

* Except the following balloons 4.0-25 mm, 4.0-30 mm and 4.0-40 whose maximum recommended pressure is 14 bars.



Balloon Diameter (mm)	Balloon Length (mm)				
	10	15	20	30	40
1.5	RAF 15-10	RAF 15-15	RAF 15-20		
2.0	RAF 20-10	RAF 20-15	RAF 20-20	RAF 20-30	RAF 20-40
2.5	RAF 25-10	RAF 25-15	RAF 25-20	RAF 25-30	RAF 25-40
3.0	RAF 30-10	RAF 30-15	RAF 30-20	RAF 30-30	RAF 30-40
3.5	RAF 35-10	RAF 35-15	RAF 35-20	RAF 35-30	RAF 35-40
4.0	RAF 40-10	RAF 40-15	RAF 40-20	RAF 40-30	RAF 40-40

Distributed by



4, Passage Saint Antoine
92508 Rueil-Malmaison Cedex - France
Tel. +33 (0)1 41 39 01 92
Fax. +33 (0)1 41 39 01 93



No contractual document - Rev01-09/2006

