

Responder AED Pro

Automated External Defibrillator

GE introduces its professional AED offering – the Responder AED Pro. With its ECG color display, manual- override feature, and multiple-rescue memory, it is the most technologically advanced device available to healthcare and rescue professionals for emergency treatment of sudden cardiac arrest and continuous cardiac monitoring.

Designed after our flagship, feature-rich, Responder AED, the Pro offers customers the same ground-breaking rescue options, programmable functions, and monitoring capabilities. Patented RescueReady® technology assures first time, every time shock delivery. Intuitive and comprehensive voice and text prompts guide the user through the rescue.

- ECG color display
- Manual-override operation
- Configurable biphasic energy protocols and detection rates
- Patented RescueReady technology includes daily, weekly, and monthly self-tests of battery, internal electronics, software and pads
- Optional continuous ECG patient monitoring
- Multiple-rescue internal memory storage
- Optional rechargeable battery
- Infrared data communications



Defibrillator	
Operation	Semi-Automatic and manual
Waveform	Biphasic truncated exponential
Energy (J) Range Available	150J - 300J, selectable via ServiceLink; 105J - 360J range of delivered energy
Number of Protocols	5
Shock Time Cycle (typical)	Approximately 55 seconds for three shock series
Control Buttons	Automated (shock) and manual
Voice Prompts	Comprehensive voice instructions guide user through rescue process
Display	Displays written instructions to guide user through rescue process, SmartGauge Battery Status Indicator, Service Indicator, Pads Indicator, Text Display, ECG Display
Visible Indicators	RescueReady Status Indicator
Audible Alerts	Voice Prompt, System Alert
Synchronized Shock	Built-in automatic synchronization feature
Pacemaker Pulse Detection	Yes
Programmable	Yes, via ServiceLink®
Pediatric Capability	Yes
Warranty	7 years

Defibrillator Electrodes	
Minimum Combined Surface Area	228 cm ²
Extended Length of Lead Wire Supplied	1.3 m
Type	Self-checking, pre-connected to the AED
Shelf Life	Adult, pre-gelled, self-adhesive, disposable defibrillation electrodes 2 years

Each Responder AED Pro package includes one (1) defibrillator, one (1) disposable IntelliSense lithium battery, one (1) pair of defibrillation electrodes, one (1) Service CD (with service manual, ServiceLink customization software, and ServiceLink instructions), one (1) RescueLink CD (with RescueLink event-review software), and one (1) operator's manual.

GE Medical Systems Information Technologies GmbH
Munzinger Strasse 3
79111 Freiburg, Germany
Tel. +49 761 4543 0 • Fax +49 761 4543 233

www.gehealthcare.com



GE imagination at work

Battery	
Type	IntelliSense® lithium
Warranty	1 years or 12 hours of use, whichever occurs first

Automated Self-Tests	
Daily	Battery, pads (presence and function), internal electronics, no energy charge, SHOCK/CONTINUE button, and software
Weekly	Battery, pads (presence and function), internal electronics, partial energy charge, SHOCK/CONTINUE button, and software
Monthly	Battery, pads (presence and function), internal electronics, full energy charge cycle, SHOCK/CONTINUE button, and software

Event Documentation	
Internal Memory	60 minutes ECG data with event annotation
Multiple Rescue	Yes
Playback	Viewable via RescueLink® software via PC
Communications	Infrared Data Transfer to PC via adapter
Clock Synchronization	Precise synchronization of AED time to dispatch time

Physical Dimensions	
Height	8,7 cm (3,4 in.)
Width	27 cm (10,6 in.)
Depth	31,8 cm (12,5 in.)
Weight	3,1 kg (6,8 lb.)

Options	
ECG cable for continuous patient monitoring	
Rechargeable IntelliSense lithium battery	

Certification	
	CE marking indicating compliance with the Medical Device Directive 93/42/EEC

©2006 General Electric Company – All rights reserved.

General Electric Company reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation.

GE, GE Monogram and Responder™ are trademarks of General Electric Company.

ServiceLink®, SmartGauge®, STAR®, IntelliSense®, RescueReady®, RescueLink® and RHYTHMx® are trademarks and registered trademarks of Cardiac Science, Inc.

GE Medical Systems Information Technologies GmbH, a General Electric company, going to market as GE Healthcare.