



HEMOCHRON® *Response*

Whole Blood Coagulation System



For all your hemostasis management needs

For more than 30 years, HEMOCHRON® has been the standard for performing whole blood coagulation monitoring at the point-of-care. HEMOCHRON Response V2.0 includes an easy-to-read screen and data management features, such as two level quality control lockout to enhance regulatory compliance plus optional dosing.





HEMOCHRON Response V2.0 with Data Manager V3.0



QUALITY CONTROL

- 2 level Liquid Quality Control* (LQC) lockout.
- 2 level Electronic System Verification* (ESV) lockout.
- Set LQC and ESV interval from 1 to 1,080 hours.
- Stores 300 QC results per well.
- Temperature Verification Tube (TVT).
 - Enter test result.
 - Save in database.

*Option for 0, 1 or 2 level QC lockout by well



HEMOCHRON RESPONSE DATA MANAGER V3.0 (SUPPLIED WITH RESPONSE V2.0)

- Allows supervisor to create and lock OID table.
 - Need PC and HRDM to modify table.
 - Updates captured in audit trail.
- Operator Identification (OID).
 - Set OID from 3 to 9 digits, or 0 for any length from 1 to 9 digits.
 - Store up to 504 OIDs.
 - Upload OID table to Response.
- Patient Identification (PID).
 - Set PID length from 3 to 9 digits or 0 for any length from 1 to 9 digits.
- Create and store up to 9 User Notes.
 - Upload to Response.
- Download results to a PC.
- View patient or QC data in Excel spread sheet.




OTHER FEATURES

- Set LQC and ESV intervals from 1 to 1,080 hours.
- Bar coded tubes.
 - Test name and expiration date automatically read when using a bar coded tube.
 - Results are stamped with date and time.
- Attach up to 2 User Notes per record.
- Results displayed as whole blood, plasma equivalent, or INR for PT assay only.
- Store 600 patient tests per well.
- Stored results can be reviewed by:
 - Test type
 - PID
 - OID
 - Date
- Internal printer
 - Select on, off or automatic printing of data.
- Languages
 - English, Spanish, French, Italian, German, Portuguese.

RXDX® ANALYSIS AND DOSING

- Optional module for heparin and protamine dosing.
- Program activated by serial number.
- Stores 10 cases and 32 tests per case.
- QC lockouts will not be enabled until a case is completed.



HEMOCHRON® and RxDx® are registered trademarks of ITC.



8 Olsen Avenue • Edison, NJ 08820 USA
 tel: 732.548.5700 • fax: 732.248.1928
 ITC Europe • Strada Rivoltana • 20090 Rodano (MI) ITALY
 tel: +39.02.9532.0196 • fax: +39.02.9532.0276
 www.itcmed.com

a subsidiary of Thoratec Corporation