

Retrieval Technology

AndraBasket & AndraBasket Micro




- The Power Of Extraction
- High Performance
- Ease Of Use



DELIVERING LESS INVASIVE FUTURE TECHNOLOGY

AndraBasket & AndraBasket Micro

1. Atraumatic tip to easy position the basket
2. High-Performance Nitinol struts ensure best capture abilities
3. Highly flexible PU shaft
4. Supportive sleeve avoids kinking, even in difficult procedures
5. Ease of Use handle device
6. Luer Port



* The AndraBasket has been specifically designed to meet the high requirements during extraction. It makes it your first choice extraction device

Order Information AndraBasket with extraction Handle:

Basket Shape	French Size	Basket Diameter	No. Struts	Sheath Length (cm)	Order Codes
Helical	2.5	12 mm	4	90	ABT-H-4-2.5-12-090H
Helical	3	15 mm	4	90	ABT-H-4-3-15-090H
Helical	4	17 mm	4	90	ABT-H-4-4-17-090H
Helical	2.5	12 mm	4	120	ABT-H-4-2-12-120H
Helical	3	15 mm	4	120	ABT-H-4-3-15-120H
Helical	4	17 mm	4	120	ABT-H-4-4-17-120H

HELICAL SHAPE



Order Information AndraBasket Micro with extraction Handle: AndraBasket Micro has a braided shaft

Basket Shape	French Size	Basket Diameter	No. Struts	Sheath Length (cm)	Order Codes
Helical	2.0	12 mm	4	90	ABT-H-4-2-12-090H

High Quality Products Made In Germany!

Your local Distributor:



Andramed GmbH
 Silvanerstraße 8
 D - 72766 Reutlingen
 Germany
www.andramed.com
 Tel.: +49-7127-972290
 Fax: +49-7127-972291

CE 0483
 Rev. 01/07