


Retrieval Technology

AndraSnare & AndraSnare Micro




- **Your First Choice Safety Device**
- **High Performance**
- **Ease of Use**



DELIVERING LESS INVASIVE FUTURE TECHNOLOGY

AndraSnare & AndraSnare Micro



1. Atraumatic Nitinol loop with gold-plated tungsten for optimal visibility
2. Pre-angled tip eases navigation
3. Radiopaque platinum inner marker ring ensures nearest positioning of introducer
4. Shapable introducer for high performance
5. Supportive sleeve strongly connected to the hub avoids kinking
6. Luer Connector
7. "SLIGHT" insertion tool for easy loading and minimizing backbleeding
8. Extra "BIGG" torque device to ease manipulation
9. One piece nitinol wire enhances the near 1:1 torque for better manipulation

* The AndraSnare & AndraSnare Micro have been specifically designed to meet the high requirements during extraction. It makes it your first choice safety device

Order Information AndraSnare:

Product Code	Loop-Diameter (mm)	Length (cm)	Introducer Sheath (cm)	Introducer Sheath (F)
AS-5	5mm (3-5mm vessel)	125	110	4
AS-10	10mm (6-10mm vessel)	125	110	4
AS-15	15mm (11-15mm vessel)	125	110	4
AS-20	20mm (16-20mm vessel)	125	110	4
AS-25	25mm (21-25mm vessel)	125	110	5
AS-30	30mm (26-30mm vessel)	125	110	5
AS-35	35mm (31-35mm vessel)	125	110	5

Product Code	Loop-Diameter (mm)	Length (cm)	Introducer Sheath (cm)	Introducer Sheath (F)
AS-25/65	25mm (21-25mm vessel)	65	50	5
AS-25/150	25mm (21-25mm vessel)	150	135	5

Order Information AndraSnare Micro made for small vessels:

Product Code	Loop-Diameter (mm)	Length (cm)	Introducer Sheath (cm)	Introducer Sheath (F)
ASM-2	2mm (1-2mm vessel)	175	150	3-2.3
ASM-4	4mm (3-4mm vessel)	175	150	3-2.3
ASM-7	7mm (5-7mm vessel)	175	150	3-2.3

High Quality Products Made In Germany!

Your local Distributor:



Andramed GmbH
 Silvanerstraße 8
 D - 72766 Reutlingen
 Germany
www.andramed.com
 Tel.: +49-7127-972290
 Fax: +49-7127-972291

CE 0483
 Rev.01/07