

COMPLIANCE TABLE

Pressure in Bars	Ø 2.5 mm	Ø 3.0 mm	Ø 3.5 mm	Ø 4.0 mm
6	2.40	2.87	3.34	3.82
7	2.45	2.93	3.42	3.90
<b>Nominal Pressure</b> → 8	<b>2.50</b>	<b>3.00</b>	<b>3.50</b>	<b>4.00</b>
9	2.55	3.05	3.55	4.07
10	2.60	3.10	3.60	4.14
✓ 11	<b>2.65</b>	<b>3.15</b>	<b>3.64</b>	<b>4.21</b>
12	2.70	3.20	3.69	4.28
13	2.74	3.25	3.74	4.36
14	2.78	3.31	3.80	4.43
15	2.82	3.35	3.85	4.50
<b>Rated Burst Pressure*</b> → 16	<b>2.86</b>	<b>3.40</b>	<b>3.90</b>	<b>4.55</b>
17	2.90	3.44	3.94	4.60
18	2.94	3.49	3.99	4.65
19	2.98	3.53	4.04	4.70
20	3.03	3.57	4.08	4.75



CORONARY STENT  
PRE-MOUNTED ON  
RAPID EXCHANGE  
PTCA POLYAMIDE  
CATHETER

- Usable length: 140 cm
- Shaft diameter:
  - Polytube 2.6 F (0.87 mm)
  - Hypotube 1.9 F (0.63 mm)
- Minimum guiding catheter 5F (0.058" / 1.47mm)
- Maximum guide wire diameter: 0.014" (0.356 mm)
- Radio-opaque markers: one at each end of stent
- Pre-mounted stent diameter: ≤ 1 mm
- Percentage of shortening: 5% at 3.0 mm
- Degree of recoil: 2.5% at 3.0 mm
- Material composition: 316L stainless steel
- Metallic surface area: 17.9% at 3.0 mm

Stent length (mm)	Balloon diameter (mm)	Reference	Deployment pressure	Rated burst pressure
7 mm	2.5 mm	LICM 25-7 S	11 bars	16 bars
	3.0 mm	LICM 30-7 S	11 bars	16 bars
	3.5 mm	LICM 35-7 S	11 bars	16 bars
10 mm	2.5 mm	LICM 25-10 S	11 bars	16 bars
	3.0 mm	LICM 30-10 S	11 bars	16 bars
	3.5 mm	LICM 35-10 S	11 bars	16 bars
13 mm	2.5 mm	LICM 25-13 S	11 bars	16 bars
	3.0 mm	LICM 30-13 S	11 bars	16 bars
	3.5 mm	LICM 35-13 S	11 bars	16 bars
	4.0 mm	LICM 40-13 S	11 bars	16 bars
16 mm	2.5 mm	LICM 25-16 S	11 bars	16 bars
	3.0 mm	LICM 30-16 S	11 bars	16 bars
	3.5 mm	LICM 35-16 S	11 bars	16 bars
	4.0 mm	LICM 40-16 S	11 bars	16 bars
19 mm	2.5 mm	LICM 25-19 S	11 bars	16 bars
	3.0 mm	LICM 30-19 S	11 bars	16 bars
	3.5 mm	LICM 35-19 S	11 bars	16 bars
	4.0 mm	LICM 40-19 S	11 bars	16 bars
22 mm	2.5 mm	LICM 25-22 S	11 bars	16 bars
	3.0 mm	LICM 30-22 S	11 bars	16 bars
	3.5 mm	LICM 35-22 S	11 bars	16 bars
	4.0 mm	LICM 40-22 S	11 bars	14 bars
28 mm	2.5 mm	LICM 25-28 S	11 bars	16 bars
	3.0 mm	LICM 30-28 S	11 bars	16 bars
	3.5 mm	LICM 35-28 S	11 bars	16 bars
	4.0 mm	LICM 40-28 S	11 bars	14 bars
38 mm	3.0 mm	LICM 30-38 S	11 bars	16 bars
	3.5 mm	LICM 35-38 S	11 bars	14 bars
	4.0 mm	LICM 40-38 S	11 bars	14 bars

No contractual document - Rev0 - 04/2005

distributed by



HEADQUARTERS

4, Passage Saint Antoine  
92508 Rueil-Malmaison Cedex - France  
Tel : (33) 01 41 39 01 92  
Fax : (33) 01 41 39 01 93

SUBSIDIARIES

HEXACATH ASIA PACIFIC	HEXACATH ESPAÑA	HEXACATH Oy
HEXACATH BENELUX	HEXACATH ITALIA	HEXACATH SUISSE
HEXACATH GmbH	HEXACATH INDIA	HEXACATH UK

